



Store at: +2+8°C.

Presentation:


Cod. SU024 CONT: R1 1 x 100 mL + R2 10 x 10 mL + CAL 1 x 5 mL.

**Procedure**

**Quantitative determination of lactate.**  
Only for *in vitro* use in clinical laboratory (IVD)

**TEST SUMMARY**

Lactate is oxidized by lactate oxidase (LO) to pyruvate and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), which under the influence of peroxidase (POD), 4-aminophenazone (4-AP) and 4-chlorophenol form a red quinone compound:



The intensity of the color formed is proportional to the lactate concentration in the sample<sup>1</sup>.

**REAGENTS COMPOSITION**

<b>R.1 (Buffer)</b>	PIPES pH 7.5 4-chlorophenol	50 mmol/L. 4 mmol/L.
<b>R.2 (Enzymes)</b>	Lactate Oxidase (LO) Peroxidase (POD) 4-Aminophenazone (4-AP)	800 U/L. 2000 U/L. 0.4 mmol/L.
<b>Lactate Cal</b>	Lactate aqueous primary calibrator	10 mg/dL.

**REAGENT PREPARATION AND STABILITY**

Working Reagent (WR): Dissolve (→) the contents of one vial R.2 (Enzymes) in 10 mL of R.1 (Buffer).

Cap and mix gently to dissolve contents. (WR) is stable 1 month at 2-8°C or 1 week at room temperature (15-25°C).

All the components of the kit are stable until the expiration date on the label when stored at 2-8°C, protected from light and contamination prevented during their use.

Do not use reagents over the expiration date.

**Lactate Cal:** Once open is stable up to 1 month when stored tightly closed at 2-8°C, protected from light and contamination prevented during their use.

**Signs of Reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm. ≥ 0.18

All the reagents of the kit are stable up to the end of the indicated month and year of expiry. Store tightly closed at 2-8°C. Do not use reagents over the expiration date.

**SPECIMEN**

Heparinized plasma. Free of hemolysis<sup>1</sup>.

Specimen must be placed on a refrigerator and separated the plasma within 15 min; because the blood cells will metabolise glucose to lactic acid. Once the plasma is separated, lactate is stable.

**MATERIAL REQUIRED BUT NOT PROVIDED**

- Spectrophotometer or colorimeter measuring at 505 nm.
- Matched cuvettes 1.0 cm. light path.

**General laboratory equipment.**

**TEST PROCEDURE**

**1. Assay Conditions**

- Wavelength : ..... 505 nm.
  - Cuvette: ..... 1 cm light path.
  - Temperature ..... 37° / 15-25°C.
- Adjust the instrument to zero with distilled water.
  - Pipette into a cuvette:

	Blank	Calibrator	Sample
WR (mL.)	1.0	1.0	1.0
Calibrator <sup>(note 1-2)</sup> (µL.)	--	25	--
Sample (µL.)	--	--	25

- Mix and incubate for 5 minutes at 37°C or 10 minutes at room temperature (15-25°C).
- Read the absorbance (A) of the samples and calibrator, against the Blank. The colour is stable at least 30 minutes.

**CALCULATIONS**

$$\text{Lactate (mg/dL.)} = \frac{(A)\text{Sample}}{(A)\text{Standard}} \times 10 \text{ (Calibrator conc.)}$$

Conversion Factor. mg/dL. x 0.111 = mmol/L.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of the procedure, Normal and Pathological. If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

**Serum controls are recommended for internal quality control. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.**

**REFERENCE VALUES**

4,5 – 19,8 mg/dL.                      0,5-2,2 mmol/L.

(These values are for orientation purpose).

**It is suggested that each laboratory establish its own reference range.**

**CLINICAL SIGNIFICANCE**

Lactate is a metabolic intermediary, originated in the lactic fermentation from glucose, which accumulates during high intensity exercise as a result of the associated increase in glycolytic activity. The formation of ATP is linked to the generation of lactate and H<sup>+</sup>. If fatigue develops, the increased levels of lactate correlate with the reduction of force<sup>1,4,5</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENT PERFORMANCE**

- **Measuring Range:**

From detection limit of 0.39 mg/dL. to linearity limit of 150 mg/dL., under the described assay conditions.

If results obtained were greater than linearity limit, dilute the sample ½ with NaCl 9 g/L. and multiply result by 2.

- **Precision:**

Mean (mg/dL)	Intra-assay n= 20		Inter-assay n= 20	
	11.1	21.8	11.4	22.1
SD	0.24	0.25	0.36	0.54
CV	2.14	1.16	3.12	2.47

- **Sensitivity:** 1 mg/dL. = 0.01A

- **Accuracy:** Results obtained GPL reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.998

Regression Equation: y=0.9979x + 1.2518

The results of the performance characteristics depend on the analyzer used.

**INTERFERING SUBSTANCES**

- Intravenous injection of epinephrine, glucose, bicarbonate, or other infusions that modify the acid-base balance, causing an elevation in lactate. Avoid using hemolyzed samples<sup>1</sup>.
- A list of drugs and other interfering substances with lactate determination has been reported by Young et. al<sup>2,3</sup>.

**NOTES**

- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

- Gau N. Lactic acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1040-1042 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.