



Store at: +2+8°C.

Presentation:

Cod. SE024 100 Test.

Procedure

Diagnostic reagent for qualitative measurement of Brucella antibodies.**Only for in vitro use in clinical laboratory (IVD)****TEST SUMMARY**

The Rose Bengal is a slide agglutination test for the qualitative and semi-quantitative detection of antibodies anti-Brucella in human and animal serum. The stained bacterial suspension agglutinates when mixed with samples containing specific IgG or IgM antibodies present in the patient sample.

REAGENTS COMPOSITION

Rose Bengal 5 mL	Brucella Abortus suspension, strain 544, in lactate buffer 1 mol/L, phenol 5 g/L, Rose Bengal, pH 3.6.
Control (+) 1 mL	Animal serum, with an antibody anti-Br. abortus concentration \geq 50 IU/mL. Sodium azide 0.95 g/L.
Control (-) 1 mL	Animal serum. Sodium azide 0.95 g/L.

PRECAUTIONS

Phenol: Toxic (T). R24/25: Toxic in contact with skin and if swallowed. R34: Causes burns. S28.2: After contact with the skin, wash immediately with plenty of water. S45: In case of accident, seek medical advice immediately.

Good laboratory safety practices should be followed when handling laboratory reagents or human samples.

REAGENT PREPARATION AND STABILITY

All the components are ready to use.

Do not use reagents over the expiration date.

Do not freeze: frozen reagents could change the functionality of the test.

If appear particles and turbidity do not use.

All the reagents of the kit are stable up to the end of the indicated month and year of expiry. Store tightly closed at 2-8°C., Do not use reagents over the expiration date.

CALIBRATION

The Rose Bengal sensitivity is calibrated against the 2^o International Preparation of anti-Brucella abortus from NIBS (UK)(WHO).

SPECIMEN

Fresh serum. Stable 8 days at 2-8°C or 3 months at -20°C.

The samples with particles or fibrin should be centrifuged to eliminate them.

Do not use haemolized or lipemic samples.

Discard contaminated specimen.

MATERIAL REQUIRED BUT NOT PROVIDED

Mechanical rotator with adjustable speed at 80-100 r.p.m.

General laboratory equipment

TEST PROCEDURE**Qualitative method**

1. Allow the reagents and sample to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50 μ L of the sample and one drop of each Positive and Negative control into separate circles on the slide test.
3. Shake the Rose Bengal reagent gently before using and add a drop of this reagent next to the sample to be tested.
4. Mix both drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
5. Rotate the slide with a mechanical rotator at 80-100 r.p.m. for 4 minutes. False positive results could appear if the test is read later than four minutes.

Semi-quantitative method

1. Make two fold dilutions of the sample in 9 gr/L Saline Solution.
2. Test each dilution as described in the qualitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates an antibody anti-Brucella concentration equal or greater than 25 IU/mL.

The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

CALCULATIONS

The approximate antibody concentration in the patient sample is calculated as follows: 25 x anti-Brucella Titer = IU/mL

QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

Serum controls Brucella are recommended for internal quality control. Each laboratory should establish its own Quality Control scheme and corrective actions.

REFERENCE VALUES

Up to 25 IU/mL.

It is suggested that each laboratory establish its own reference range.

CLINICAL SIGNIFICANCE

Brucella diagnostic may be assessed either by micro organism isolation in blood or stools, or by titration of specific antibodies in the patient serum. The reagent, because of its formulation in an acid buffer, is reactive with both IgG and IgM antibodies and very useful for the diagnosis of chronic individuals which present a high level of IgG antibody, difficult to be detected by the reference tube method (Wright).

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENT PERFORMANCE

- Analytical sensitivity:
25 (\pm 5) IU/mL, under the described assay conditions
- Prozone effect:
No prozone effect was detected up to 1000 IU/mL.
- Diagnostic sensitivity: 93.3 %
- Diagnostic specificity: 100 %

INTERFERING SUBSTANCES

Interferences:

- Hemoglobin (10 g/L), lipemia (10 g/L), rheumatoid factors (300 IU/mL) do not interfere. Bilirubin interferes at 2.5 mg/dL.
- Other substances may interfere⁵.

BIBLIOGRAPHY

1. Young E J. Clinical Infectious Diseases 1995; 21: 283-290.
2. Alton GC. Techniques for Brucellosis Laboratory INRA Paris, 1988.
3. Ariza J. Current Opinion in Infectious Diseases 1996; 9: 126-131.
4. Comité mixto FAO/OMS de expertos en Brucellosis. WLD Health Org Tech Rep Ser 1958; 148: 1-60.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

